




ALPHAVAC

MULTIPURPOSE MECHANICAL ASPIRATION



CONTROLLED REMOVAL. MAXIMUM VERSATILITY. FLOW RESTORED.

The AlphaVac F18⁸⁵ System provides procedural time efficiency while maintaining safety and clinical efficacy.



**INDICATED FOR THE TREATMENT OF
PULMONARY EMBOLISM**

The AlphaVac F18⁸⁵ System was designed to **optimize navigation** between the pulmonary arteries.



0% MORTALITY
RATE¹

35.5%
MEAN REDUCTION IN CLOT
BURDEN, VIA MODIFIED
MILLER INDEX SCORE¹



Experience the
AlphaVac System
Alpha-Vac.com



FEATURES

| | |
|-----------------------|-------------------------|
| CANNULA | F18⁸⁵ |
| US PART NUMBER | H787253020 |
| CE PART NUMBER | H78725302CE0 |
| ACCESS | 22F |
| CANNULA SIZE | 18F |
| LENGTH | 105cm |
| FUNNEL SIZE | 11mm |
| ANGLE | 85° |

REFERENCES

1. Data on file.

ALPHAVAC MMA F18⁸⁵ SYSTEM RISK INFORMATION

US Indication for Use:

The Cannula is indicated for: the non-surgical removal of thrombi or emboli from the vasculature; aspiration of contrast media and other fluids from the vasculature.

The Cannula is intended for use in the venous system and for the treatment of pulmonary embolism.

The Handle is indicated as a vacuum source for the AlphaVac Multipurpose Mechanical Aspiration System.

CE Indications for Use:

The Cannula is indicated for: the non-surgical removal of thrombi or emboli from the pulmonary arteries; aspiration of contrast media and other fluids from the pulmonary arteries.

The Cannula is intended for the treatment of pulmonary embolism.

The Handle is indicated as a vacuum source for the AlphaVac Multipurpose Mechanical Aspiration System.

Canada Indications for Use:

The Cannula is indicated for: the non-surgical removal of thrombi or emboli from the venous vasculature; aspiration of contrast media and other fluids from the venous vasculature.

The Cannula is intended for use in the venous system.

The Handle is indicated as a vacuum source for the AlphaVac Multipurpose Mechanical Aspiration System.

Refer to Instructions for Use and/or User Manual provided with the product for complete Instructions, Warnings, Precautions, Potential Complications and Contraindications prior to use of the product.

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.

The AlphaVac F18⁸⁵ System is indicated for the treatment of pulmonary embolism (PE) within the US and EU. The AlphaVac F22²⁵ and F22¹⁸⁰ Systems are not indicated for the treatment of PE, and none of the AlphaVac Systems are indicated for the treatment of PE in Canada.

The AlphaVac System is not available in all territories. Please consult your AngioDynamics representative for availability.



AngioDynamics, the AngioDynamics logo, AlphaVac and the AlphaVac logo are trademarks and/or registered trademarks of AngioDynamics, Inc., an affiliate or subsidiary. ©2024 AngioDynamics, Inc. GL/VI/BR/2840 Rev 01 08/2024



PULMONARY ARTERIES

INTERNAL
JUGULAR

ACCESS SITES

FEMORAL VEIN

